

## FOOT-UP®

Foot-Up is a lightweight AFO that offers dynamic support for drop foot or similar complaints. The orthosis provides visible improvement in gait by supporting the foot the moment it is raised.



### INDICATIONS FOR USE

- Drop foot
- Support of dorsiflexion

### FEATURES AND BENEFITS

- Streamlined, lightweight and invisible when worn
- Breathable material, padded and comfortable to wear
- Easy to fit
- Available in black and beige
- Attached to the shoe, or available as shoeless model

### REQUIRED MEASUREMENTS

- Basic version: Circumference of ankle, measured around the malleoli
- Shoeless version: Circumference of foot immediately behind the ball of the foot

### ORDERING INFORMATION

#### FOOT-UP

Size	Ankle Circumference	Beige	Black
Medium	18 - 21cm	07811 -1	07810-1
Large	22 - 26cm	07811 -2	07810-2
X-Large	27 - 33cm	07811 -3	07810-3

#### SHOELESS VERSION\*

Size	Foot Circumference	Beige	Black
Small	< 20cm	27811-1	27810-1
Medium	20 - 24cm	27811-2	27810-2
Large	24 - 28cm	27811-3	27810-3
X-Large	> 28cm	27811-4	27810-4

\*Shoeless version item comprises of foot wrap only and must be ordered with the Foot-Up unless the patient already has one.

### ACCESSORIES

Description	Beige	Black
Foot-Up plastic insert	17811	17810

Description	Beige	Black
Foot-Up with Shoeless add-on Bundle	07811xy	07810xy

x = size of main foot-up part  
 y = size of shoeless add-on  
 e.g. 0781012 = size 1 Foot-up with a size 2 Shoeless add-on in Black