

REBOUND FOOT-UP®

Rebound Foot-Up is a lightweight ankle/foot orthosis (AFO) that offers dynamic support for drop-foot or related disorders that require dorsiflexion support. The Rebound Foot-Up supports the ankle dorsiflexion during swing phase. Dynamic support is provided by an elastic strap connected to the shoe by either a transparent plastic inlay placed between shoelaces, or by a nut and bolt fastened through the eyelets. The positioning of the fastening allows for individualising the support by providing the option of positioning the dorsiflexion support on the medial or lateral side. When using the plastic inlay, the dorsiflexion support becomes neutral.

The Rebound Foot-Up Ankle Cuff has a smart clamp that holds the cuff in place and facilitates one-handed donning and doffing.

Optional: Rebound Foot-Up Foot Wrap can be used when barefoot or in open shoes or sandals. The adjustable elastic strap is attached to the Foot Wrap using a nut and bolt fastened through the eyelets.



INDICATIONS FOR USE

- Drop-foot
- Disorders that require dorsiflexion support

FEATURES AND BENEFITS

- Dynamic and discreet support for drop-foot or related disorders that require dorsiflexion support
- Improved function, ease of use, comfort, and quality
- Allows for voluntary plantar flexion
- The elastic strap is adjustable and easily replaceable
- The Rebound Foot-Up Ankle Cuff and Rebound Foot-Up Foot Wrap are made from breathable material that can be worn comfortably for long periods of time
- The Rebound Foot-Up is hand-washable
- Positioned medial or lateral support
- Contoured to the shape of leg
- Adjustable strap length
- Smart clamp facilitates one-handed donning / doffing
- Optional Rebound Foot-Up Foot Wrap can be used without shoes or in open shoes or sandals

REQUIRED MEASUREMENTS

- Circumference of ankle, measured around the malleoli



ORDERING INFORMATION

ANKLE CUFF

Description	Size	Ankle Circumference	Part #
Rebound Foot-Up Ankle Cuff	S/M	18 - 25cm (7 - 10")	RFU50001
Rebound Foot-Up Ankle Cuff	L/LX	26 - 33cm (10 - 13")	RFU50002

FOOT WRAP

Description	Size	Foot Circumference	Part #
Rebound Foot-Up Foot Wrap	One Size	20 - 28cm (8 - 11")	RFU50101

ACCESSORIES

Description	Part #
Rebound Foot-Up Attachment Kit (includes plastic inlay, elastic strap, and 1 x nut and screw)	RFU50201
Rebound Foot-Up Replacement Nuts And Bolts (5 nuts included: 2 short and 3 long screws)	RFU50202

Description	Size	Foot Circumference	Part #
Rebound Foot-Up Bundle	S/M	18 - 25cm	RFU500011
Rebound Foot-Up Bundle	L/LX	26 - 33cm	RFU500021