

# PHILADELPHIA® TRACHEOTOMY COLLAR

The Philadelphia Tracheotomy Collar is a hospital and rehabilitation collar that is designed to maintain neutral alignment of the cervical spine after injury to the neck. The collar contains a large tracheal opening, allows carotid pulse monitoring and emergency tracheotomies.

### INDICATIONS FOR USE

- C-spine precaution for trauma patients
- Immobilisation for pre c-spine surgery
- Other conditions requiring gross immobilisation of the mid-cervical spine

### FEATURES AND BENEFITS

- Anterior and posterior rigid plastic reinforcement supports limit movement
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy
- Easily customised for individual requirements
- Does not contain natural rubber latex, non-toxic and hypoallergenic Plastazote® foam reduces skin irritation



### REQUIRED MEASUREMENTS

- Neck circumference and height from tip of chin to sternal notch

## ORDERING INFORMATION

Size	Neck circumference	6cm height	8cm height	11cm height	13cm height
X-Small	< 25cm	PHP-T2XS			
Small	25 - 33cm	PHP-T2S	PHP-T3S	PHP-T4S	PHP-T5S
Medium	33 - 41cm	PHP-T2M	PHP-T3M	PHP-T4M	PHP-T5M
Large	41 - 48cm	PHP-T2L	PHP-T3L	PHP-T4L	PHP-T5L
X-Large	> 48cm	PHP-T2XL	PHP-T3XL	PHP-T4XL	PHP-T5XL

Size	Neck circumference	Part#
Infant	15 - 20cm	PHP-TI
Paediatric	20 - 28cm	PHP-TP

## ACCESSORIES

### PHILADELPHIA COOLMAX COTTON LINER

Size	Neck circumference	Part#
Child	20 - 28cm	PHP-LCP
Small	25 - 33cm	PHP-LCS
Medium	33 - 41cm	PHP-LCM
Large	41 - 48cm	PHP-LCL

### PHILADELPHIA TRACHEOTOMY ATTACHMENTS

Size	Stabiliser	Halo
Paediatric	PHP-OSP	PHP-OHP
Adult	PHP-OS	PHP-OH